


Mining Event-related Brain Dynamics I



Scott Makeig





Institute for Neural Computation
University of California San Diego

33rd EEGLAB Workshop
Aspet, France
July 4, 2023



William Blake


Swartz Center for Computational Neuroscience, UCSD












EEGLAB History

- 1993 – ERSP (Makeig)
- 1995 – Infomax ICA for EEG (Makeig, Bell, Jung, Sejnowski)
- **1997 - EEG/ICA Toolbox (cnl.salk.edu), ITC & ERC**
- 1999 - ERP-image plotting (Jung & Makeig)
- **2000 –GUI design & name ‘EEGLAB’ (Delorme)**
- **2002 – 1st EEGLAB (sccn.ucsd.edu)**
- **2004 - 1st EEGLAB support from U.S. NIH and reference paper (Delorme & Makeig, 2004)**
- 2006 - 1st EEGLAB plug-ins, STUDY structure, and component clustering tools
- 2009+ – New toolboxes: NFT, SIFT, BCILAB, MPT, ... (Akalin Acar, Mullen, Kothe, ...)
- **2011 – EEGLAB, the most widely used EEG research environment (Henke & Halchenko)**
- 2013 – Lab Streaming Layer (LSL) (Kothe) for Mobile Brain/Body Imaging (MoBI) (Makeig)
- 2013 – *HeadIT.org* online, HED neuroinformatic tools (Bigdely-Shamlo)
- 2017 – LIMO / GLM integrated (Pernet) -- and 24rd- 26th EEGLAB Workshops ...
- **2018 – The Open EEGLAB Portal via the Neuroscience Gateway (nsgportal.org).**
- 2020 – EEGLAB 2019, BIDS integration, ICLabel, get_chanlocs, ...
- 2022 -- HED tools, NEMAR,



EEGLAB Site Visits (in 24 hours)





List of data import extensions

| Plug-in name ↕ | Version ↕ | Short plug-in description ↕ | Link ↕ | Contact ↕ | Comments ↕ |
|---|-----------|--|--|---|-------------------------------|
| MFFimport  | 1.00 | Import MFF files from the EGI company | Download  | S. Chennu  | User comments |
| ANTeepimport  | 1.10 | Import ANT .cnt data and trigger files | Download  | M. van de Velde  | User comments |
| BCI2000import  | 0.36 | Import BCI2000 data files | Download  | C. Boulay  | User comments |
| BDFimport | 1.10 | Import BDF data files | Download  | A. Delorme  | User comments |
| biopac | 1.00 | Import BIOPAC data files | Download  | A. Delorme  | User comments |
| ctfimport | 1.04 | Import CTF (MEG) data files | Download  | D. Weber  | User comments |
| erpssimport | 1.01 | Import ERPS data files | Download  | A. Delorme  | User comments |
| INSTEPascimport | 1.00 | Import INSTEP ASCII data files | Download  | A. Delorme  | User comments |
| neuroimaging4d | 1.00 | Import Neuroimaging4d data files | Download  | C. Wienbruch  | User comments |
| ProcomInfinity | 1.00 | Import Procom Infinity data files | Download  | A. Delorme  | User comments |
| WearableSensing | 1.09 | Import Wearable Sensing files | Download  | S. Millen  | User comments |
| NihonKoden | 0.10 | Import Nihon Koden M00 files (beta) | Download  | M. Miyakoshi  | User comments |
| xdfimport | 1.12 | Import files in XDF format | Download  | C. Kothe  | User comments |
| bva-io  | 1.5.12 | Import Brain Vision Analyser data files | Download  | A. Widmann  | User comments |
| Fileio  | Daily | Import multiple data files formats | Download  | R. Oostenveld  | User comments |
| Biosig  | 2.88 | Import multiple data files formats | Download  | A. Schloegl  | User comments |
| Cogniscan  | 1.1 | Import Cogniscan data files | Download  | P. Sajda  | User comments |
| NeurOne  | 1.0.3.2 | Import NeurOne data files | Download  | Support  | User comments |
| loadhdf5 | 1.0 | Load hdf5 files recorded with g.recorder | Download  | Simon L. Kappel  | User comments |

EEGLAB
EXTENSION
MANAGER

EEG as Functional Brain Imaging

Hemodynamic imaging

= imaging local brain

Energy

Direct 3-D inverse model,
but quite **slow** & **indirect**
as well as **expensive**,
very heavy & **non-portable**.

1993 -

Electromagnetic imaging

= imaging

local cortical field synchrony

3-D imaging needs head model,
but a quite **fast** & **direct** measure
of *one aspect* of cortical activity –
local spatial field coherence.


1926 -

Functional Brain Imaging using EEG

- EEG imaging is noninvasive → little ethical concern
- EEG imaging can be tolerated by most subjects
- EEG imaging has fine time resolution
- EEG imaging is lightweight / mobile / wearable
- EEG imaging is inexpensive → scalable
- EEG source imaging requires a *good* forward-problem electrical head model and inverse localization method.
- Historically, much inertia in EEG methods development

Three Modern Eras of EEG Research

We are here



Loo, Lenartowicz & Makeig, 2015


Figure 1. Relative number of PubMed citations retrieved by 'All Fields' search terms: 'EEG,' 'ERP,' and 'Brain Oscillations.' The percent of citations for each search term relative to the total number of citations returned by a search for any of the three terms is plotted relative to the other two search terms. For visual clarity, 'Brain Oscillations' citations are graphed with a green dotted line according to the Y-axis labels on the right; 'EEG' with a blue solid line and 'ERP' with a red dashed line according to the Y-axis labels on the left.

What is EEG?

- Brain electrical activity
- A small portion of *cortical* brain electrical activity
- An even smaller portion of *total* brain electrical activity

- **But *which* portion?**
- **Triggered and modulated *how*?**
- **With *what* functional significance?**

The generation and modulation of local field potentials is COMPLEX and not well studied



**Information flow is
bi-directional !**

**Local
Extracellular
Fields**

**What is an
EEG
"source" ?**

**Brain dynamics are
inherently multi-scale**

EEG (scalp surface fields)

**ECOG (larger cortical
surface fields)**


At each spatial recording scale, the signal is produced by **active partial coherence** of distributed activities at the next smaller spatial scale.

**Intracellular and
peri-cellular fields**

**Synaptic and
other trans-
membrane
potentials**

Phase cones (Freeman)






S. Anderson


Simple patterns

Complex patterns


Delta band
(1-4 Hz)
in
anesth.
animals




“Synchrony was associated with high delta-band amplitude (averaged across the recording array), whereas complex waves were associated with low average delta-band amplitude. ...




Spike rates were highest near the position and time of spirals and saddles and lowest in the presence of synchrony.”





1-cm ECoG
chan spacing




Tchoe et al., *Science Translat. Med.*, 2022

1-mm ECoG
chan spacing

Naïve 2-D interpretation of EEG signals?




Cortical EEG signal projection patterns as point processes



Cortical source current volume conduction patterns


The very broad EEG point-spread function



Single simulated parietal source →

Very broad projected scalp potentials


The very broad EEG point-spread function



Single simulated parietal source →

Very broad projected scalp potentials


The very broad EEG point-spread function



Single simulated parietal source →

Very broad projected scalp potentials


The very broad EEG point-spread function




Single simulated parietal source →

Very broad projected scalp potentials

The 'receptive field' of a bipolar EEG channel





Scalp EEG channel



Its cortical 'receptive field'

The 'receptive field' of a bipolar EEG channel!


Scalp EEG channel $A \rightarrow B$




Its cortical 'receptive field'

At time t ,

$$V_{A \rightarrow B}^{(t)} = \sum_{r \text{ in Cortex}} v_r \times W_{AB}(r)$$




Bipolar channel




Average reference channel




Phase cones (Freeman)






An Effective EEG Source






Here, what are the cortical 'effective sources'?



Scalp projection


Scalp epiphenomena !

Phenomena





**Two spatially stationary
cortical effective sources**

Epiphenomena



**Summed
scalp projection**




A close-up photograph of a human eye with a light blue iris and dark eyelashes. The eye is looking slightly to the right. Overlaid on the eye in a bold, blue, sans-serif font is the text "Who am I?".

Who
am I?

Embodied Agency

Brain processes
have evolved and function
to optimize the outcomes
of the **willed** behavior
the brain organizes
in response to
perceived & felt
challenges and opportunities.

Brains meet the challenge
of the moment
– *every moment!*



Three Aspects of Human Consciousness

Knowing - I perceive, remember, believe

Feeling - I feel, experience as feeling

Willing - I act, aim, intend

“[Humans] have *full consciousness* of the [physical] world
in **all the aspects of knowing, feeling and willing.**”

Meher Baba

EEG can be used to learn and monitor
how the brain and nervous system
support human consciousness
in ALL its aspects --

Knowing

Feeling

Willing

Electromagnetic source localization using realistic head models



NFT


Solve for forward problem using realistic head model (EM)

T1-weighted



MRI


Segmentation



Simple Map

Signal Processing

Source Estimate




EEG/MEG

Blind EEG Source Separation by Independent Component Analysis



Tony Bell,
developer of
Infomax ICA

ICA can find distinct EEG source activities -- and their 'simple' scalp maps!



Independent Component Analysis of Electroencephalographic Data

Scott Makeig
Naval Health Research Center
P.O. Box 85122
San Diego CA 92186-5122
scott@epi.lanog.nhrc.navy.mil


Anthony J. Bell
Computational Neurobiology Lab
The Salk Institute, P.O. Box 85800
San Diego, CA 92186-5800
tonyb@salk.edu

Tzyy-Ping Jung
Naval Health Research Center and
Computational Neurobiology Lab
The Salk Institute, P.O. Box 85800
San Diego, CA 92186-5800
jung@salk.edu

Tiziana J. Sejnowski
Howard Hughes Medical Institute and
Computational Neurobiology Lab
The Salk Institute, P.O. Box 85800
San Diego, CA 92186-5800
terry@salk.edu

Abstract

Because of the distance between the skull and brain and their different sensitivities, electroencephalographic (EEG) data collected from any point on the human scalp includes activity generated within a large brain area. This spatial smearing of EEG data by volume conduction does not involve significant time delays, however, suggesting that the Independent Component Analysis (ICA) algorithm of Bell and Sejnowski [1] is suitable for performing blind source separation on EEG data. The ICA algorithm separates the problem of source identification from that of source localization. First results of applying the ICA algorithm to EEG and event-related potential (ERP) data collected during a sustained auditory detection task show: (1) ICA training is insensitive to different random seeds; (2) ICA may be used to segregate obvious artifactual ERP components (eye and muscle noise, eye movements) from other sources; (3) ICA is capable of isolating overlapping ERP phenomena, including alpha and theta bursts and spatially-separable ERP components, to separate ICA channels; (4) Nonstationarities in EEG and behavioral state can be tracked using ICA via changes in the amount of residual correlation between ICA-filtered output channels.




Are EEG 'effective source' signals independent?

Independent

Domains

of Local Synchrony



The EEG Inverse Problem is Twofold

Effective source


Identification → Localization

ICA gives a model-based response to the first question:


- ***What* are the effective sources? (identification)**

And it greatly helps answer the second question:


- ***Where* do these sources originate? (localization)**



ICA separates *non-brain* effective source processes








ICA finds non-brain independent component (IC) processes ...




... separates them from the remainder of the data ...

... including IC EMG sources



... and *also* separates cortical brain IC processes



ICLabel: A crowd-sourced AI independent component classifier

L. Pion-Tonachini et al.

ICLabel

NeuroImage 198 (2019) 181–197

Newest IC

Back 1

Log Out


User: [SCCN] Luca Pion-Tonachini
Labels: 5054
Don't use the browser back button.

Leaderboard


Tutorial

Profile




011360_006



Component Time Series



IC 6 of 63
Data Var. Accounted For:
4.99%



Brain

Muscle

Eye

Heart

Next

Line Noise

Chan Noise

Other


?

Fig. 1. An IC labeling example from the ICLabel website (<https://iclabel.ucsd.edu/tutorial>), which also gives a detailed description of the features shown above. Label contributors are shown the illustrated IC measures and must decide which IC category or categories best apply. They mark their decision by clicking on the blue buttons below, and have the option of selecting multiple categories in the case that they cannot decide on one or believe the IC contains an additive mixture of sources. There is also a “?” button that they can use to indicate low confidence in the submitted label.

labeling.ucsd.edu/tutorial


ICLabel: A crowd-sourced AI independent component classifier

036669_017




ICLabel

Component Time Series




IC 17 of 68
Data Var. Accounted For:
0.47%

Epoched Data



1 dipole


2 dipole



RV: 3.7%


RV: 3.1%
DMR: 6.3


Activity Power Spectrum



labeling.ucsd.edu/tutorial

ICA in practice





Klaus Gramann
3-D Tunnel Task


A Passive Spatial Navigation Paradigm




A Passive Spatial Navigation Paradigm






‘Turner’ and ‘Nonturner’ subjects use different spatial orienting styles



Two parietal component clusters





Medial prefrontal component cluster




Turn visible

Exit visible



Clusters distinguishing Turners & Nonturners



The Beginning

EEG brain imaging etc. . . .

