


DIPFIT and model co-registration

1. Co-register electrodes with model
2. Fit components







Edit channel info -- pop_chanedit()



Channel information ("field_name"):

Channel label ("label")	HEOG	Opt. head center
Polar angle ("theta")	-42	Rotate axis
Polar radius ("radius")	0.65556	Transform axes
Cartesian X ("X")	55.7734	
Cartesian Y ("Y")	50.2186	XYZ -> polar & sph.
Cartesian Z ("Z")	-39.9051	Sph. -> polar & xyz
Spherical horiz. angle ("sph_theta")	42	Polar -> sph. & xyz
Spherical azimuth angle ("sph_phi")	-28	
Spherical radius ("sph_radius")	85	Set head radius
Channel type		Set channel types
Reference		Set reference
Index in backup 'urchanlocs' structure	68	
Channel in data array (set=yes)	<input checked="" type="checkbox"/>	

Buttons: Delete chan, Insert chan, Channel number (of 72): 68, Append chan

Plot 2-D Plot radius (0.2-1, []=auto) Nose along +X Plot 3-D (xyz)





EEGLAB development head

File Edit **Tools** Plot Study Datasets Help

- #1: EE
- Change sampling rate
- Filter the data
- Re-reference
- Interpolate electrodes
- Reject continuous data by eye
- Extract epochs
- Remove baseline
- Run ICA
- Remove components
- Automatic channel rejection
- Automatic continuous rejection
- Automatic epoch rejection
- Reject data epochs
- Reject data using ICA
- ICLabel
- Clean continuous data using ASR
- Locate dipoles using DIPFIT**

coregister()

File Edit View Insert Tools Desktop Window Help

Labels off
Electrodes
Labels on
Electrodes
Mesh off

Move right (mm)	0	Pitch (rad)	0	Resize (x)	99.05	Align fiducials
Move front (mm)	0	Roll (rad)	0	Resize (y)	99.05	Warp montage
Move up (mm)	0	Yaw (rad)	0	Resize (z)	99.05	Cancel Ok



coregister()

File Edit View Insert Tools Desktop Window Help


Labels off
Electrodes
Labels on
Electrodes
Mesh off

Move right (mm)	0	Pitch (rad)	0	Resize (x)	99.05	Align fiducials
Move front (mm)	0	Roll (rad)	0	Resize (y)	99.05	Warp montage
Move up (mm)	0	Yaw (rad)	-1.571	Resize (z)	99.05	Cancel Ok

Head model and settings

- Component dipole coarse fit
- Component dipole fine fit
- Component dipole plot
- Component dipole autofit
- Distributed source Leadfield matrix
- Distributed source component modelling
- Source reconstruction of ERP

EEG.dipfit structure



```
>> EEG.dipfit
```

```
ans =
```

```
hdmfile: [1x76 char]
```

```
mrifile: [1x71 char]
```

```
chanfile: [1x83 char]
```

```
chansel: [1x33 double]
```

```
coordformat: 'spherical'
```

```
model: [1x33 struct]
```

```
current: 32
```

```
vol: [1x1 struct]
```

```
leadfield: []
```

```
coord_transform: [0 0 -1.570796 100 76 90.87264 1 1 1]
```




Rotation

Translation

Scaling





Plot scalp maps in 3D




Go through co-registration
in the same way as
with dipfit co-registration






2D scalp map for IC 12



3D scalp map for IC 12







DIPFIT and model co-registration


1. Co-register electrodes with model
2. Fit components





Patch of Cortex Acting as a Dipole





Fit equivalent dipoles

EEGLAB development head


File Edit **Tools** Plot Study Datasets Help

#1: EE




- Change sampling rate
- Filter the data
- Re-reference
- Interpolate electrodes
- Reject continuous data by eye
- Extract epochs
- Remove baseline
- Run ICA
- Remove components
- Automatic channel rejection
- Automatic continuous rejection
- Automatic epoch rejection
- Reject data epochs
- Reject data using ICA
- ICLabel
- Clean continuous data using ASR
- Locate dipoles using DIPFIT

Head model and settings


- Component dipole coarse fit
- Component dipole fine fit
- Component dipole plot
- Component dipole autofit**
- Distributed source Leadfield matrix
- Distributed source component modelling
- Source reconstruction of ERP



Computing residual variance



$$r = \frac{\sum(x_i - \tilde{x}_i)^2}{\sum x_i^2}$$



Delorme A, Palmer J, Onton J, Oostenveld R, Makeig S. Independent EEG sources are dipolar. PLoS One. 2012;7(2):e30135. doi: 10.1371/journal.pone.0030135. Epub 2012 Feb 15. PMID: 22355308; PMCID: PMC3280242.

Scroll through dipoles





Figure 3

File Edit View Insert Tools Desktop Window

1 dipoles:


- Plot one
- Keep|Next
- Next
- Prev
- Keep|Prev
- 1
- Comp: 3
- RV: 2.91%
- X tal: 4
- Y tal: -81
- Z tal: 33
- cuneus R**
- Display:**
- Mesh on
- Tight view
- Sagittal view
- Coronal view
- Top view
- No controls

Figure 3



Tools Desktop Window Help

- Next
- Prev
- Keep|Prev
- 1
- Comp: 3
- RV: 2.91%
- X tal: 4
- Y tal: -81
- Z tal: 33
- cuneus R
- Display:**
- Mesh on
- Tight view
- Sagittal view
- Coronal view
- Top view
- No controls


Desikan-Killiany Atlas




Right Cuneus







68 brain areas






Visualizing ICA component clusters







Distributed source localization in DIPFIT





- Surface source model: Colin27 (with Desikan-Kilianny atlas)
- Surface source model: Use Brainstorm ICBM152 (with Desikan-Kilianny atlas)
- Volumetric source model: LORETA-KEY
- Volumetric source model: AFNI with TTatlas+tirc atlas (Fieldtrip)
- Custom source model






- Head model and settings
- Component dipole coarse fit
- Component dipole fine fit
- Component dipole plot
- Component dipole autofit
- Distributed source Leadfield matrix**
- Distributed source component modelling
- Source reconstruction of ERP









- Surface source model: Colin27 (with Desikan-Kilianny atlas)
- Surface source model: Use Brainstorm ICBM152 (with Desikan-Kilianny atlas)
- ✓ Volumetric source model: LORETA-KEY
- Volumetric source model: AFNI with TTatlas+tlrc atlas (Fieldtrip)
- Custom source model






Distributed source localization (eloreta or LCMV beamforming)



- Head model and settings
- Component dipole coarse fit
- Component dipole fine fit
- Component dipole plot
- Component dipole autofit
- Distributed source Leadfield matrix
- Distributed source component modelling**
- Source reconstruction of ERP



Exercise



- Load ***eeglab_data_epochs_ica.set*** from the ***sample_data*** folder of EEGLAB
- Remove EOG channels
- In Dipfit menu, select BEM model
- Coarse fit all dipoles
- Fine fit dipole number 5
- Compute Leadfield matrix
- Compute distributed source for component 3

